

Digital Multi-Position Magnetic Stirrer series

MS3-MaxiStir4P

4 Stirring point quantity

MS3-MaxiStir8P

8 Stirring point quantity

MS3-MaxiStir12P

12 Stirring point quantity

MS3-MaxiStir8PT

8 Stirring point quantity
with temperature control

Multi-position Magnetic Stirrer is designed for the high—efficiency experiments and has 4/8/12 stirring bar configurations to choose. It can use laboratory space efficiently. DC brushless motor makes the model, low speed but stable, high speed but powerful with low noise. Stainless steel is coated by silicone resin film disk. Non-slip and anti-corrosion.

Digital display stirring speed at real-time, stirring bars operate synchronously and provide uniform agitation reaction conditions.

Stainless steel is coated by silicone resin film pad. Non-slip and anti-corrosion.

Model MS3-MaxiStir8PT applies special heating technology, surface temperature up to 120°C



MS3-MaxiStir4P



MS3-MaxiStir8P



MS3-MaxiStir12P



MS3-MaxiStir8PT

Specification	MS3-MaxiStir4P	MS3-MaxiStir8P	MS3-MaxiStir12P	MS3-MaxiStir8PT
Stirring Point Quantity	4	8	12	8
Size of Stirring Bars	φ8 x 46mm	φ8 x 46mm	φ8 x 46mm	φ8 x 46mm
Max. Capacity	400ml x 4	400ml x 8	400ml x 12	400ml x 8
Digital Display	Speed Display	Speed Display	Speed Display	Temperature and speed
Speed Range	200-1200rpm	± 200-1200rpm	± 200-1200rpm	200-1200rpm
Temp. Control Range	/	/	/	R.T. +5°C-120°C
Platform Material	Stainless steel with silicone pad	Stainless steel with silicone pad	Stainless steel with silicone pad	Stainless steel with silicone pad
Motor	DC brushless motor	DC brushless motor	DC brushless motor	DC brushless motor
Voltage	AC 220V, 50/60Hz	AC 220V, 50/60Hz	AC 220V, 50/60Hz	AC 220V, 50/60Hz
Power	20W	20W	20W	400W
Fuse	250V, 1A, φ5x20	250V, 1A, φ5x20	250V, 1A, φ5x20	250V, 5A/3A, φ5x20
Dimension	W.115 x D.480 x H.50mm	W.205 x D.480 x H.50mm	W.295 x D.480 x H.50mm	W.205 x D.480 x H.50mm
Net weight	3kgs	4kgs	5kgs	5.3kgs