

# MS3-MaxiRocker3D

## 3D Gyraryory Rocker

3D Gyraryory Rocker is an economical 3D gyraryory rocker with perfect 3-dimensional motion for gentle but thorough mixing. It is an ideal choice for mixing small sample. It suits for blood collection tubes and centrifuge tubes, applied for all blot related applications such as gel staining and destaining. The concave pad can firmly keep the tubes.

Gentle but thorough mixing.

Suits for blood collection tubes and centrifuge tubes.

Fixed speed and tilt angle.

Safe for use in cold rooms and incubators

MS3-MaxiRocker3D is available for 1.5ml x 60, 3ml x 30, 7 ml x 21, 15ml x 17.

Specification	MS3-MaxiRocker3D
Speed Range	20rpm (fixed)
Shaking Orbit	3D / 20° (fixed)
Max. Load	0.8 kgs
Voltage	AC 220V / AC 110V, 50/60Hz
Power	45W
Dimension	W.240 x D.170 x H.150mm
Net Weight	1.6kgs

