

MS3-MaxiRocker80 / Rocking shaker **MS3-MaxiRocker80D / Rocking shaker Double platform**

MS3-MaxiRocker80 and MS3-MaxiRocker80D Rocking shaker is a desktop oscillator that can be used directly in the greenhouse or in the incubator. It is widely used in the fixation of various electrophoresis gels, cocmassie blue staining, ooscillating shaking during decolorization, silver nitrate staining fix, dye, develop and other laboratory experiments

Platform is made of stainless steel and natural rubber, which is acid and alkali resistant, non-deformation, easy-celan and solution spill-prevention.


Real time display of speed and diminishing time at the control panel

Simple and beautiful appearance.

MS3-MaxiRocker80 and MS3-MaxiRocker80D Products includes flat platform with non-slip rubber mat.



MS3-MaxiRocker80



MS3-MaxiRocker80D

Specification	MS3-MaxiRocker80	MS3-MaxiRocker80D
Shaker Speed Range	5-80 rpm	5-80 rpm
Shaker Mechanism	25mm (up and down) / 10°	25mm (up and down) / 10°
Digital Timer Control Range	1min-99h59min	1min-99h59min
Max. Loading capacity	5kgs	5kgs
Platform Dimension	L. 280 x W.228 mm	L. 280 x W.228 mm Double platform
Voltage	AC 220-230V, 50/60Hz	AC 220-230V, 50/60Hz
Power	50W	50W
Fuse	250V, 1A , Ø5x20	250V, 1A , Ø5x20
Dimensions, WxDxH	284x264x134 mm	284x264x248 mm
Net weight	0,5,2 kgs	6 kgs